

A Complete Pipetting Solution

Axygen® Axypet® Pro Single- and Multi-channel Pipettors and Axxygen MultiRack Pipet Tip System



Axygen Axypet Pro pipettors and the Axxygen MultiRack tip system have been engineered to provide the highest levels of comfort, accuracy, and precision.

The lightweight construction, contoured handle, and 4-digit counter of Axxygen Axypet Pro pipettors ensure comfortable pipetting. All feature smooth plunger movement and extremely low pipetting forces to reduce wrist strain and fatigue. Pushbuttons on the pipettors are color-coded to match the racks of Axxygen MultiRack tips to enable easy size identification and reduce the risk of cross-contamination.

The Axxygen MultiRack tip system is an interchangeable, complete system offered in bulk, racked, filter, and reload packaging styles to meet the variety of needs in the lab setting. The MultiRack tip system offers standard and extended-length tips allowing volumes from 0.1 μL to 1250 μL . For racked and reload packaging styles, the tips are pre-loaded into easy-load, color-coded deck inserts. With a small benchtop footprint and reusable rack, this system minimizes environmental impact and lab space. The stackable design saves critical bench space. Axxygen MultiRack tip system products are made in the USA and meet ISO standards. All tips are free of detectable RNase, DNase, DNA, and pyrogens.

Axygen Axypet Pro Pipettors

- Ergonomic handle design fits both the right and left hand
- Convenient one-hand volume setting with auto-lock system prevents accidental volume change during work
- 4-digit counter provides greater precision and is positioned for visibility at all times
- Single-channel pipettors are compatible with very narrow tubes
- Full autoclavability and UV-resistance helps prevent sample contamination
- The multi-channel shafts retract individually for easier tip loading and ejection
- Unique ejector lever system enables effortless tip ejection
- Easy in-lab calibration helps you quickly adjust your pipettor to various liquid types
- Individually tested and supplied with a certificate of quality

Axygen MultiRack Pipet Tips

- Reusable rack
- Small footprint and stackable design for maximum bench space
- Extended tips offered in 10, 200, 300, and 1000 μL sizes (1000 μL extended length allows up to 1250 μL volume)

Compatibility Chart



Axygen Axypet Pro Pipettor		Axygen MultiRack Pipet Tips									
		Non-filtered				Filtered					
Cat. No.	Description	10 µL	200 µL	300 µL	1000 µL	10 µL	20 µL	100 µL	200 µL	300 µL	1000 µL
AP-2-P	0.1-2 µL single-channel pipettor	■				■					
AP-10-P	0.5-10 µL single-channel pipettor	■				■					
AP-20-P	2-20 µL single-channel pipettor		■				■				
AP-50-P	5-50 µL single-channel pipettor		■					■			
AP-100-P	10-100 µL single-channel pipettor		■					■			
AP-200-P	20-200 µL single-channel pipettor		■						■		
AP-1000-P	100-1,000 µL single-channel pipettor				■						■
AP-8-10-P	0.5-10 µL, 8-channel pipettor	■				■					
AP-8-50-P	5-50 µL, 8-channel pipettor		■								
AP-8-200-P	10-200 µL, 8-channel pipettor		■	■					■	■	
AP-8-300-P	30-300 µL, 8-channel pipettor			■						■	
AP-12-10-P	0.5-10 µL, 12-channel pipettor	■				■					
AP-12-50-P	5-50 µL, 12-channel pipettor		■								
AP-12-200-P	10-200 µL, 12-channel pipettor		■	■					■	■	
AP-12-300-P	30-300 µL, 12-channel pipettor			■						■	

Ordering Information

Products may not be available in all markets.

Axygen® Axypet® Pro Pipettors

Single-channel

Cat. No.	Volume Range (µL)	Accuracy (%)	Precision (%)	Non-filtered Tip (µL)	Qty/Cs
AP-2-P	0.1 - 2	±40.0 to ±1.5	≤12.0 to ≤0.7	10	1
AP-10-P	0.5 - 10	±4.0 to ±1.0	≤2.8 to ≤0.4	10	1
AP-20-P	2 - 20	±3.5 to ±0.8	≤1.5 to ≤0.3	200	1
AP-50-P	5 - 50	±3.5 to ±0.8	≤1.3 to ≤0.3	200	1
AP-100-P	10 - 100	±3.0 to ±0.8	≤1.0 to ≤0.2	200	1
AP-200-P	20 - 200	±2.0 to ±0.6	≤0.7 to ≤0.2	200	1
AP-1000-P	100 - 1000	±2.5 to ±0.6	≤0.6 to ≤0.2	1000	1

8-channel

AP-8-10-P	0.5 - 10	±10.0 to ±2.0	≤8.0 to ≤1.2	10	1
AP-8-50-P	5 - 50	±4.0 to ±1.6	≤2.5 to ≤0.6	200	1
AP-8-200-P	10 - 200	±5.0 to ±1.0	≤3.0 to ≤0.6	200/300	1
AP-8-300-P	30 - 300	±3.0 to ±1.0	≤1.5 to ≤0.6	300	1

12-channel

AP-12-10-P	0.5 - 10	±10.0 to ±2.0	≤8.0 to ≤1.2	10	1
AP-12-50-P	5 - 50	±4.0 to ±1.6	≤2.5 to ≤0.6	200	1
AP-12-200-P	10 - 200	±5.0 to ±1.0	≤3.0 to ≤0.6	200/300	1
AP-12-300-P	30 - 300	±3.0 to ±1.0	≤1.5 to ≤0.6	300	1


Axygen Axypet Pro Starter Kit

Cat. No.	Description	Qty/Cs
AP-STR-KIT-P	Axygen Axypet Pro Starter Kit Includes: 4 single-channel pipettors: 0.5 to 10 µL (Cat. No. AP-10-P), 2 to 20 µL (Cat. No. AP-20-P), 20 to 200 µL (Cat. No. AP-200-P), 100 to 1000 µL (Cat. No. AP-1000-P); universal linear stand for four single-channel pipettors (Cat. No. AP-USD4-P); Axypet MultiRack pipet tips: 10 µL (Cat. No. MR-10XT-R), 200 µL (Cat. No. MR-200-R), 1000 µL (Cat. No. MR-1000XT-R)	1

Accessories

AP-USD	Pipettor stand for one single-channel, 8- or 12-channel pipettor	1
AP-USD4-P	Universal linear stand, red, for four pipettors (four single-channel, two 8- or 12-channel pipettors)	1
AP-USD4	Universal linear stand, transparent, for four pipettors (four single-channel, two 8- or 12-channel pipettors)	1
AP-USD6	Universal linear stand, transparent, for six pipettors (six single-channel, three 8- or 12-channel pipettors)	1





Axygen® MultiRack Pipet Tips

Bulk

Cat. No.	Description	Quantity
MR-10XT	10 µL MultiRack extended length pipet tip, bulk	1,000/pk, 10 pks/cs
MR-200	200 µL MultiRack pipet tip, bulk	1,000/pk, 10 pks/cs
MR-1000XT	1000 µL MultiRack extended length pipet tip, bulk	768/pk, 10 pks/cs

Racked

MR-10XT-R	10 µL MultiRack extended length pipet tip, racked	960/pk, 5 pks/cs
MR-200-R	200 µL MultiRack pipet tip, racked	960/pk, 5 pks/cs
MR-300-R	300 µL MultiRack pipet tip, racked	960/pk, 5 pks/cs
MR-300XT-R	300 µL MultiRack pipet tip, extended length, racked	768/pk, 5 pks/cs
MR-1000XT-R	1000 µL MultiRack extended length pipet tip, racked	768/pk, 5 pks/cs

Racked, Sterile

MR-10XT-R-S	10 µL MultiRack extended length pipet tip, racked, sterile	960/pk, 5 pks/cs
MR-200-R-S	200 µL MultiRack pipet tip, racked, sterile	960/pk, 5 pks/cs
MR-300-R-S	300 µL MultiRack pipet tip, racked, sterile	960/pk, 5 pks/cs
MR-300XT-R-S	300 µL MultiRack extended length pipet tip, racked, sterile	768/pk, 5 pks/cs
MR-1000XT-R-S	1000 µL MultiRack extended length pipet tip, racked, sterile	768/pk, 5 pks/cs

Axygen MultiRack Reload System

Cat. No.	Description	Quantity
MR-10XT-RL	10 µL MultiRack extended length pipet tip, reload	960/pk, 4 pks/cs
MR-200-RL	200 µL MultiRack pipet tip, reload	960/pk, 4 pks/cs
MR-300-RL	300 µL MultiRack pipet tip, reload	960/pk, 4 pks/cs
MR-1000XT-RL	1000 µL MultiRack extended length pipet tip, reload	960/pk, 4 pks/cs

Axygen MultiRack Filter Pipet Tips

Racked, Filter, Sterile

Cat. No.	Description	Quantity
MRF-10XT-R-S	10 µL MultiRack extended length filtered pipet tip, racked, sterile	960/pk, 5 pks/cs
MRF-20-R-S	20 µL MultiRack filtered pipet tip, racked, sterile	960/pk, 5 pks/cs
MRF-100-R-S	100 µL MultiRack filtered pipet tip, racked, sterile	960/pk, 5 pks/cs
MRF-200NX-R-S	200 µL NX MultiRack filtered pipet tip, racked, sterile	960/pk, 5 pks/cs
MRF-200XT-R-S	200 µL MultiRack extended length filtered pipet tip, racked, sterile	768/pk, 5 pks/cs
MRF-300XT-R-S	300 µL MultiRack extended length filtered pipet tip, racked, sterile	768/pk, 5 pks/cs
MRF-1000XT-R-S	1000 µL MultiRack extended length filtered pipet tip, racked, sterile	768/pk, 5 pks/cs

Axygen MultiRack Pipet Tips with Axxygen Maxymum Recovery® Surface

Racked, Sterile, Axxygen Maxymum Recovery Surface

Cat. No.	Description	Quantity
MR-10XT-L-R-S	10 µL MultiRack extended length pipet tip, Maxymum Recovery surface, racked, sterile	960/pk, 5 pks/cs
MR-200-L-R-S	200 µL MultiRack natural pipet tip, Maxymum Recovery surface, racked, sterile	960/pk, 5 pks/cs
MR-300-L-R-S	300 µL MultiRack pipet tip, Maxymum Recovery surface, racked, sterile	960/pk, 5 pks/cs
MR-300XT-L-R-S	300 µL MultiRack extended length pipet tip, Maxymum Recovery surface, racked, sterile	768/pk, 5 pks/cs
MR-1000XT-L-R-S	1000 µL MultiRack extended length pipet tip, Maxymum Recovery surface, racked, sterile	768/pk, 5 pks/cs



Axygen® MultiRack Pipet Tips with Axygen Maxymum Recovery® Surface (Continued)

Racked, Filter, Sterile, Axygen Maxymum Recovery Surface

Cat. No.	Description	Quantity
MRF-10XT-L-R-S	10 µL MultiRack extended length filtered pipet tip, Maxymum Recovery surface, racked, sterile	960/pk, 5 pks/cs
MRF-20-L-R-S	20 µL MultiRack filtered pipet tip, Maxymum Recovery surface, racked, sterile	960/pk, 5 pks/cs
MRF-100-L-R-S	100 µL MultiRack filtered pipet tip, Maxymum Recovery surface, racked, sterile	960/pk, 5 pks/cs
MRF-200NX-L-R-S	200 µL NX MultiRack filtered pipet tip, Maxymum Recovery surface, racked, sterile	960/pk, 5 pks/cs
MRF-200XT-L-R-S	200 µL MultiRack extended length filtered pipet tip, Maxymum Recovery surface, racked, sterile	768/pk, 5 pks/cs
MRF-300XT-L-R-S	300 µL MultiRack extended length filtered pipet tip, Maxymum Recovery surface, racked, sterile	768/pk, 5 pks/cs
MRF-1000XT-L-R-S	1000 µL MultiRack extended length, filtered pipet tip, Maxymum Recovery surface, racked, sterile	768/pk, 5 pks/cs

Axygen MultiRack Tip Station

Cat. No.	Description	Quantity
MR-10XT-TS	10 µL MultiRack extended length pipet tip station	192/layer, 5 layers/stack, 960/stack, 5 stacks/cs
MR-10XT-TS-S	10 µL MultiRack extended length pipet tip station, sterile	192/layer, 5 layers/stack, 960/stack, 5 stacks/cs
MR-200-TS	200 µL MultiRack pipet tip station	192/layer, 5 layers/stack, 960/stack, 5 stacks/cs
MR-200-TS-S	200 µL MultiRack pipet tip station, sterile	192/layer, 5 layers/stack, 960/stack, 5 stacks/cs

Accessories

Cat. No.	Description	Quantity
MR-R-A	Empty rack for Axygen MultiRack 10 µL extended length, 200 µL, 200 µL NX, 300 µL	5/pk, 1 pk/cs
MR-R-B	Empty rack for Axygen MultiRack 200 µL extended length, 300 µL extended length, and 1000 µL extended length	5/pk, 1 pk/cs



Warranty/Disclaimer: Unless otherwise specified, all products are for research use or general laboratory use only.* Not intended for use in diagnostic or therapeutic procedures. Not for use in humans. These products are not intended to mitigate the presence of microorganisms on surfaces or in the environment, where such organisms can be deleterious to humans or the environment. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications. *For a listing of US medical devices, regulatory classifications or specific information on claims, visit www.corning.com/resources.

Corning's products are not specifically designed and tested for diagnostic testing. Many Corning products, though not specific for diagnostic testing, can be used in the workflow and preparation of the test at the customers discretion. Customers may use these products to support their claims. We cannot make any claims or statements that our products are approved for diagnostic testing either directly or indirectly. The customer is responsible for any testing, validation, and/or regulatory submissions that may be required to support the safety and efficacy of their intended application.

CORNING

Corning Incorporated
Life Sciences

www.corning.com/lifesciences

NORTH AMERICA

t 800.492.1110
t 978.442.2200

ASIA/PACIFIC

Australia/New Zealand
t 61 427286832

Chinese Mainland
t 86 21 3338 4338

India

t 91 124 4604000

Japan

t 81 3-3586 1996

Korea

t 82 2-796-9500

Singapore

t 65 6572-9740

Taiwan

t 886 2-2716-0338

EUROPE

CSEurope@corning.com

France

t 0800 916 882

Germany

t 0800 101 1153

The Netherlands

t 020 655 79 28

United Kingdom

t 0800 376 8660

All Other European Countries

t +31 (0) 206 59 60 51

LATIN AMERICA

grupoLA@corning.com

Brazil

t 55 (11) 3089-7400

Mexico

t (52-81) 8158-8400